


OCLUSORES

DUCTUS PICCOLO




El Amplatzer Piccolo™ Occluder, que ofrece un cierre comprobado de PDA para pacientes de 700 g y más, brinda la fuerza para ocluir conductos pequeños, al tiempo que minimiza la protrusión en la arteria pulmonar circundante y la aorta descendente.

- Los discos de retención de perfil bajo y el tornillo final están diseñados para minimizar la protrusión hacia la aorta y la arteria pulmonar.
- Diseño simétrico para abordaje pulmonar o aórtico Malla de tejido apretado para minimizar la derivación residual inmediatamente después de la colocación.
- Cable de entrega flexible para facilitar la implementación.
- Capacidad de ser recapturado y redistribuido fácilmente para una colocación óptima.

DUCTUS PICCOLO

Available Device Sizes

9 Sizes Available



NOTE: Length [C] is measured from inside proximal disc to inside distal disc.

Model / Reorder Number Sizing is noted as: 9-PDAP-(waist)-(length)-L	A Waist Diameter (mm)	C Waist Length (mm)	B Disc Diameter (mm)
9-PDAP-03-02-L	3	2	4.00
9-PDAP-03-04-L	3	4	4.00
9-PDAP-03-06-L	3	6	4.00
9-PDAP-04-02-L	4	2	5.25
9-PDAP-04-04-L	4	4	5.25
9-PDAP-04-06-L	4	6	5.25
9-PDAP-05-02-L	5	2	6.50
9-PDAP-05-04-L	5	4	6.50
9-PDAP-05-06-L	5	6	6.50



Do not implant if:

- PDA length is less than 3 mm
- PDA diameter is more than 4 mm at the narrowest portion